

Anigen Rapid *New product*

FeliD-3 Ag

General Information

Purpose : Detection of Feline Parvo, Corona and Giardia antigen

Principle : Immunochromatographic assay

Specimen : Feline feces

Materials provided :

- Anigen Rapid FPV/FCoV/Giardia test device
- Assay diluents tube for FPV/FCoV Ag (P)
- Assay diluents tube for Giardia Ag (G)
- Disposable swab
- Disposable dropper




Special Features :

- Able to detect the pathogen in cats
- One-step testing procedure : Quick & Accurate
- High sensitivity and specificity
 - Sensitivity : FPV Ag 97 % vs. RT-PCR, FCoV Ag 95.5 % vs. PCR
 - Specificity : FPV Ag 98.5 % vs. RT-PCR, FCoV Ag 98.8% vs. PCR



Test Procedure

[Proper sample amount collected by a swab]



Ordering Information

Cat. No.	Description	Type	Pack size
RC1208DC	FeliD-3 Ag	Quattro Device	5 Tests/Kit